


# QUARTZ GUIDED INFRARED HEATERS



## Thermo Riko's

Quartz guided infrared heaters provide the ideal method for heating your sample in any application that requires clean, high temperature, high speed and heating in gas atmosphere by means of infrared wavelength irradiation. This patented technology was developed 30 years ago and more than 600 units have been installed in Japan.

## High Temperature

There are 3 different power selections for an Infrared lamp, 1kW, 2kW and 3kW. 1kW will be able to heat up a sample typically up to 1,300°C. 2kW will be able to heat a sample up to 1,500°C. Detailed specifications are described in the following pages.

## High Speed

Specially designed infrared lamp emits a near Infrared region of wavelengths (mainly 900nm) and the wavelength is focused through the gold coated spheroidal reflector lens. The infrared light is focused on the surface of the upper end of the Quartz rod and then the Quartz rod transmit the wavelength with minimum loss of the total reflection. Thus, very high

density infrared radiation will be emitted from the other end of the rod. The maximum heating speed can be 150°C/sec or more depending on the model.

## No Contact

Your sample surface will be heated up by Infrared radiation which comes through the Quartz rod without loss of radiation intensity. The sample will be located right at the front of the tip of the Quartz rod and will not require direct contact to the tip of the rod, because the heater won't heat up your sample by means of thermal contact.

## No outgas


From the heat source. The Infrared heater's heat source is a light bulb which is located far from the sample. In the case when your sample is under vacuum pressure, outgassing from the heat source (typically filament or sample holder) contaminates your analysis. The Infrared heater provides the cleanest method for your sample heating under vacuum environment. Furthermore, as this is non-contact, turning off the lamp will not provide any residual heat around your sample. Thus, the cool-

ing down rate is much faster than other heaters. From the max. temperature to room temperature can be within minutes, not hours.

## Oxygen environment

The Quartz and the vacuum seal are the only materials which are exposed to the vacuum. This allows experiments or processes under any gas environments. For example, heat up samples in Oxygen environment is possible.

Quartz Guided Infrared ray



# UHV

Clean heating in the UHV applications is essential, especially for analytical applications like the XPS (X-ray photoelectron spectrometer). The GVH series uses a Sapphire viewing port on a standard CF35 flange for Ultrahigh Vacuum sealing.

Applications

- Clean heating of a sample in UHV in XPS or SEM.
- Non destructive sample cleaning for surface analysis

UHV



1400°C

No contact

No outgas

Ultimate temperature : 1,400°C  
UHV Compatible  
Ultimate pressure: 5 x 10<sup>-11</sup> mbar

Principle of heater, GVH




# GVH

# GVL

# HV

| Model Name           | GVH198                                               | GVH298  | GV198                                               | GVL298  | GVL398  |
|----------------------|------------------------------------------------------|---------|-----------------------------------------------------|---------|---------|
| Lamp power           | 1kW                                                  | 2kW     | 1kW                                                 | 2kW     | 3kW     |
| Ultimate temperature | 1,200°C                                              | 1,400°C | 1,300°C                                             | 1,500°C | 1,600°C |
| Max. Heating rate    | 1 °C/sec                                             |         | 100 ~ 150°C/sec                                     |         |         |
| Heating size         | ø 20 mm (std), (min. ø5. max. ø28mm)                 |         |                                                     |         |         |
| Vacuum seals         | Copper Gasket (viewing port)                         |         | Viton O-ring                                        |         |         |
| Leak rate            | 1.33 x 10 <sup>-10</sup> Pa•m³/s or less             |         | 1.33 x 10 <sup>-8</sup> Pa•m³/s or less             |         |         |
| Ultimate pressure    | 5 x 10 <sup>-11</sup> mbar (5 x 10 <sup>-9</sup> Pa) |         | 5 x 10 <sup>-9</sup> mbar (5 x 10 <sup>-7</sup> Pa) |         |         |
| Cooling              | Air cooling fan(s) and Water cooling                 |         |                                                     |         |         |
| Water flow rate,     | 1 l/min                                              | 2 l/min | 1 l/min                                             | 2 l/min |         |
| Mounting flange      | CF35 (ø70mm OD) Conflat© flange                      |         |                                                     |         |         |



The GVL series provides a high speed heating rate. The maximum temperature can be 1,500°C with the GVL298, 2kW model and 1,500°C, can be achieved within a minute.

Applications

- Ultra high speed heating of Si, SiC and Graphene
- Formation of films or oxides crystals in oxygen atmosphere
- Laser ablation and temperature desorption spectroscopy, PLD system and rapid annealing system...etc.

High speed

1500°C


No contact

No outgas

Ultimate temperature : 1,500°C  
Heating rate: max. 150°C/sec  
Ultimate pressure: 5 x 10<sup>-9</sup> mbar



Principle of heater, GV/GVL



## Accessories


### Programmable Temperature Controller (incl. in a standard package)



#### Features

- AC100V 50/60Hz main
- Program 16 patterns/16 segments
- Max. temperature : 1,700°C
- Sensor: works with R-type T/C
- Analog out: 1 x 0~5V
- Size: 480 x 149 x 320mm

### R-Type thermocouple with linear drive, RS250V (optional)



#### Features

- Extremely thin sheathed:  $\phi 0.5$
- R type thermocouple ( $\sim 1,600^{\circ}\text{C}$ )
- Linear travel: 25mm @0.2 res.
- Mounting flange: CF16
- Leak rate:  $1.33 \times 10^{-10} \text{ Pa}\cdot\text{m}^3/\text{sec}$
- K type available ( $\sim 1,200^{\circ}\text{C}$ )


### Cooling water flow switch, DFS3R (included in a standard package)



#### Features

- Monitor cooling water flow
- Interlock signal to controller
- Flow rate: 0.5 ~ 3 L/s adj.


### Infrared thermal sensor, IR2S (optional)



#### Features

- No contact monitoring
- Range: 400 ~ 2,000°C
- Detector: IrGaAs
- Working distance: 200mm (std)

## GA Series - Atmospheric applications




Quartz guided infrared heaters are available for applications under atmospheric pressure - GA Series.

| Model name        | <b>GA298</b>   | <b>GA198</b> |
|-------------------|----------------|--------------|
| Lamp power        | 2kW            | 1kW          |
| Ultimate temp.    | 1,400°C        | 1,200°C      |
| Max. heating rate | 150°C/sec.     | 100°C/sec.   |
| Heating Area      | 20mm dia. max. |              |
| Water flow rate   | 2 L/s          | 1 L/s        |

## Examples of installation of the heater



Vertical (or inverted) mount  
with a standard  
SD series linear drive



Mount in any position with an optional HD100 linear drive

Contact:

**THERMO RIKO CO., LTD**

Mitaka High-Tech Center, 8-7-3 Shimorenjaku  
Mitaka-city, Tokyo 181-0013, JAPAN  
website: <http://www.kagaku.co.jp/thermo/>  
email: [sekigai@thermo-r.co.jp](mailto:sekigai@thermo-r.co.jp)

Sales Contact (Europe):

**Japan Vacuum Instruments**

30, Rundle Road, Liverpool, United Kingdom  
website: <http://thermo-riko.co.uk/>, <http://jpvacinst.co.uk/>  
email: [info@jpvacinst.co.uk](mailto:info@jpvacinst.co.uk)